

HIGH LEVEL OF 11-6102 -5202 PROFESSIONAL 05-750 SERVICES AND 50-750-750 SUPPORT



Why do clients choose us?

Arvand swiss is a leading manufacturer of infection prevention products Our aim is to provide high level of professional sale services and support that adds value and relevance to your staff and needs; we will deliver ongoing valued consultancy to prepare for the future growth of health care field.

Arvand Swiss is a customer centered company who provide a professional, friendly and individual service to every customer focusing on their needs.



Swi\$eal

High performance thermosealer with an innovative design

- Roll holder
- Welding band of 12mm
- According to ISO11607-2
- Acoustic and visual signal for sealing
- Easiness of paper entree
- Safety system to avoid burning of material



Swi∮eal					
Technical Features					
Code	Dimension(cm)	Power(watt)	Max valtage(Volt)	Seal Width(mm)	Seal Length(mm)
330501006	47*30*30	110	220	12	300

SWISS CUT

- Ergonomic design
- Useful length: 50, 70, 90 (cm)

SWISS CUT					
Technical Features					
Code	Product Description	Cut Width(mm)			
330599001	Arvand cut 50	450			
330599002	Arvand cut 70	550			
330501003	Arvand cut 90	650			

AR 8000 SEALCUT

- Innovative machine to produce medical pouches
- Conforms to ISO 11607-2
- With a reel dispenser and heat sealing device
- Equipped with LCD, USB port for traceability
- Maximum flexibility
- Cost effective
- Compact dimension
- Programmable

AR 8000 SEALCUT				
Technical Features				
Code	Product Description	Dimension(cm)		
4113002	AR 8000	740*355*240		

Automatic Sealing Machine

- Automatic Sealing Machine Conforms to ISO11607- 2 - With an integrated printer

Automatic Sealing Machine		
CODE	PRODUCT DESCRIPTION	
4113001	Automatic Sealing Machine	








RESISTOMETER

- Design based on EN ISO 18472 requirements
- Ready to use program for each type of test based on ISO 11140 for true positive and true negetive test
- Allow to define program with different temeperature and time requested by customer
- Capable to perform consecutive tests cycle after initial warm up
- Perform tests needed in type 1 to type 6 chemical indicators

- Equipped with different Alarms and warnings to inform user of any non conformity with standard requirements

- Print test result

RESISTOMETER				
Technical Features				
Code	Product Description	Dimension(cm)		
100575	RESISTOMETER	130*120*85		



MEDICAL PACKAGING VALIDATION LABORATORY

Arvand Swiss Medical Packaging Validation Laboratory is complying with the quality system regulation and ISO 11607 requirements. Arvand swiss comply requirements for the competence of testing and calibration according to ISO/IEC 17025.

we do the material, process and equipment qualification which includes:

- Accelerated Aging
- Barrier properties
- Seal & Tensile Strength test
- Installation qualification, Operational qualification and performance qualification
- Mechanical properties

AUTOCLAVE LABORATORY

Arvand Swiss autoclave Test and Calibration Laboratory has the most advanced technology, most sophisticated measuring equipment and employs qualified personnel specializing at calibration and quality control services in the field of medical autoclaves.

We are honored that our teamwork is based on ISO 9001 and Arvand swiss comply requirements for the competence of testing and calibration according to ISO/IEC 17025.

The scope of the Arvand Swiss Calibration Laboratory includes:

- Table-top Autoclave
- Hospital Steam sterilization
- Laboratory autoclave
- pharmaceutical Autoclave



We do quality control of Autoclave according to EN 13060, EN 285, ISO 17665, and IEC 61010



